	T			T	T
	(EN) Test cassette	Blood volume and test time	New: Performance characteristics according to	If peripheral circulation is impaired, capillary blood	
Falzgitter	EN TEST CASSELLE	The meter requires approximately 0.3 µL of blood (1 µL	EN ISO 15197:2013	might not be a true reflection of the physiological	Consult package insert
420x297 / 70x37	Read this package insert and the User's Manual of the	(microlitre) = 1 thousandth of a millilitre) per blood	System accuracy:	blood glucose level. This may apply in the following	A 0
Falzschema C	Accu-Chek Mobile blood glucose meter before testing	glucose test. The test takes approximately 5 seconds	System accuracy results for glucose concentration less	circumstances: Severe dehydration as a result of diabetic ketoacidosis or due to a hyperglycaemic	related notes in the package
	your blood glucose with this test cassette. The User's	(measuring time depends on the blood glucose concentration).	than 100 mg/dL (less than 5.55 mmol/L)	hyperosmolar non-ketotic syndrome, hypotension,	insert accompanying this
	Manual contains all the information you need to		within within within ± 5 mg/dL ± 10 mg/dL ± 15 mg/dL	shock, decompensated heart failure NYHA class IV or	product.
	perform a test. If you have any questions, contact your customer support and service centre.	Storing and using test cassettes properly	(within (within (within	peripheral arterial occlusive disease. • You may use blood with a haematocrit of 25 to 55 %.	Temperature limitation (store at)
	Castomer Support and Souther South	Test cassettes which are not stored or used	\pm 0.28 mmol/L) \pm 0.56 mmol/L) \pm 0.83 mmol/L)		at)
	This package insert features the following 2 symbols:	properly can lead to incorrect test results. Incorrect test results can cause the wrong therapy	120/156 154/156 155/156	Reagent composition	Use by
	↑ This symbol indicates a possible risk of injury or	recommendation to be made and thus produce	(76.9 %) (98.7 %) (99.4 %)	Minimum content per cm ² at time of manufacture	Use-by period of the test
	of damage to your health.	serious adverse health effects.	System accuracy results for glucose concentration	Mutant variant of quinoprotein glucose dehydrogenase (Mut. Q-GDH 2, modified	cassette after opening the
	This symbol draws your attention to important information.	Always store the test cassettes in the unopened	equal to or greater than 100 mg/dL (equal to or greater than 5.55 mmol/L)	variant of EC 1.1.5.2), acinetobacter 6.7 U	foil-sealed plastic container:
	information.	plastic container.	within ± 5 % within ± 10 % within ± 15 %	spec.	90 days
		Store the test cassettes at a temperature between +2 and +30 °C in a dry place away from direct	216/444 410/444 442/444	Pyrroloquinoline quinone 0.40 μg	All components of the pack
		sunlight.	(48.6 %) (92.3 %) (99.5 %)	Bis-(2-hydroxyethyl)-(4-hydroximino- cyclohexa-2,5-dienylidene)-ammonium 8.5 μg	domestic waste. Discard
	Intended use	If you store the test cassette in a refrigerator, leave	System accuracy results for glucose concentration	chloride	used test cassettes
	The test cassette is intended for quantitatively	the unopened plastic container to stand at an ambient temperature. Only remove the test cassette	between 29.3 mg/dL (1.6 mmol/L) and 497.0 mg/dL	2,18-phosphomolybdic acid, sodium salt 88 μg	according to local regulations.
	measuring blood glucose in fresh capillary blood.	once the plastic container has warmed up to ambient	(27.6 mmol/L)	Stabiliser 0.17 mg	Manufacturer
	The test cassette may only be used with Accu-Chek	temperature. This prevents condensation from	within ± 15 mg/dL or within ± 15 %	Non-reactive ingredients 2.1 mg	Wallulacturei
	Mobile meters and may only be used outside of the body.	forming in the test cassette. The test areas are sensitive to humidity. Only	(within \pm 0.83 mmol/L or within \pm 15 %)	Last update	REF Catalogue number
	 The system comprised of meter and test cassette is 	transport the test cassettes in the unopened plastic	597/600 (99.5 %)	2015-01	III
	only suitable for self-testing. People with diabetes	container or in the meter. When you open the plastic	Repeatability:	Customer Support and Service Centre	LOT Batch code
TITEL (EN) (RU)	can use this system to self-test their blood glucose.	container, you must use up the test cassette within 90 days (use-by period). After a period of more than	Mean value	Tastomer Support and Service Centre	In vitro diagnostic medical
	 The system must not be used to diagnose or rule out diabetes. 	90 days, the test cassette can deliver incorrect test	[mg/dL] 43 75 128 214 326	客戶服務熱線:+852-2485 7512(辦公時間)	In vitro diagnostic medical device
ACCU-CHEK®		results.	[mmol/L] 2.4 4.2 7.1 11.9 18.1 Standard deviation	www.accu-chek.com.hk	GTIN Global Trade Item Number
	Additional information	 Do not use test cassettes if the plastic container or foil cover is damaged. 	[mg/dL] 2.0 2.3 3.1 4.3 7.1	Hong Kong	UTIN
Mobile	Self-testing is not a substitute for visits to your healthcare professional. You must receive proper	 If you take a partly used test cassette out of the 	[mmol/L] 0.1 0.1 0.2 0.2 0.4	Enquiry hotline: +852-2485 7512 (office hours) www.accu-chek.com.hk	C € 0088 This product fulfils the
07141254	instruction from a qualified healthcare professional	meter, keep it in a dry place away from light.	Coefficient of variation		requirements of the European Directive 98/79/EC on in vitro
	before you start self-testing your blood glucose.	When you perform a test, the temperature must be	[%] — — 2.4 2.0 2.2	References	diagnostic medical devices.
TITEL	Your healthcare professional will determine the appropriate blood glucose range jointly with you.	between +10 and +40 °C.	Intermediate precision:	 D'Orazio et al.: Approved IFCC Recommendation on Reporting Results for Blood Glucose (Abbreviated); 	
EN RU	Keep the blood glucose monitoring system with	Test principle	Mean value	Clinical Chemistry 51:9, 1573–1576, 2005	
	all its components out of reach of children under 3	Each test area contains reagents. When blood is applied to the test area, the glucose dehydrogenase enzyme	[mg/dL] 44 126 390	[2] American Diabetes Association: Standards of	
AKKY-YEK®	years. There is a risk of suffocation if small parts (e.g. covers, caps or similar objects) are	(Mut. Q-GDH 2, EC 1.1.5.2) reacts with the blood glucose.	[mmol/L] 2.4 7.0 21.6	Medical Care in Diabetes–2012. <i>Diabetes Care</i> 35 (Suppl. 1), S11–S63, 2012	
Мобайл	swallowed.	The subsequent chemical reaction changes the colour of	Standard deviation	[3] IDF Clinical Guidelines Task Force. Global guideline	
WIOGaux	All components of the pack can be discarded in	the test area. The meter registers this colour change and converts it into a blood glucose value.	[mg/dL] 2.0 2.8 7.8	for Type 2 diabetes. Brussels: International	
07141254	domestic waste. Discard used test cassettes		[mmol/L] 0.1 0.2 0.4 Coefficient of variation	Diabetes Federation, 2005 [4] Definition and diagnosis of diabetes mellitus	
	according to local regulations. This test cassette delivers results that	Control test Carry out regular tests with the Accu-Chek Mobile	[%] 4.6 2.2 2.0	and intermediate hyperglycemia: report of a	
	correspond to blood glucose concentrations in	control solutions to ensure that your meter and test		WHO/IDF consultation. WHO, Geneva 2006	
	plasma as per the recommendation of the	cassette are functioning properly and that you are	Performance assessment by the user: A study evaluating glucose values from fingertip	(ISBN 92 4 159493 4, ISBN 978 92 4 159493 6)	
	International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) [1]. Therefore, your	carrying out tests correctly. Read the package insert for the control solutions and the User's Manual for your	capillary blood samples obtained by 104 lay persons		
	meter displays blood glucose values that refer to	meter.	showed the following results:		
	plasma although you always apply whole blood to	Performance characteristics of the	 For glucose concentrations less than 100 mg/dL (less than 5.55 mmol/L), 100 % of the test results 		
	the test area.	Accu-Chek Mobile system	were within \pm 15 mg/dL (within \pm 0.83 mmol/L) of		
	The normal fasting glucose level for a non-diabetic	The Accu-Chek Mobile system complies with the	the results obtained through laboratory testing.		
	adult is below 100 mg/dL (5.6 mmol/L). A criterion for the diagnosis of diabetes in adults is a fasting glucose	requirements of EN ISO 15197:2013 (In vitro diagnostic- test systems – Requirements for blood glucose	 For glucose concentrations equal to or greater than 100 mg/dL (equal to or greater than 5.55 mmol/L), 		
	level of 126 mg/dL or higher (7.0 mmol/L or higher)	monitoring systems for self-testing in managing	96.8 % of the test results were within \pm 15 % of the		
	confirmed in two tests [2, 3, 4]. Adults with a fasting	diabetes mellitus).	results obtained through laboratory testing.		
	glucose level between 100 and 125 mg/dL (5.6 and 6.9 mmol/L) are defined as having impaired fasting	Calibration and traceability: The system (meter and	Sources of error which may produce incorrect		
	glucose (prediabetes) [2]. Other diagnostic criteria for	test cassette) is calibrated with whole blood containing various glucose concentrations as a calibrator. The	test results		
	diabetes exist. Consult your healthcare professional to	reference values are obtained using the hexokinase	Incorrect test results can cause the wrong therapy recommendation to be made and thus produce		C € 0088
	determine if you have diabetes or not.	method which is calibrated using the ID-GCMS method. The ID-GCMS method as the method of highest	serious adverse health effects.		IN VITRO DIAGNOSTIC
	Contents of the pack	metrological quality (order) is traceable to a primary	Parenteral administration of galactose and		MEDICAL DEVICE
	1 or 2 test cassettes	NIST standard. Using this traceability chain, the test	galactosemia can lead to falsely elevated test results.		© 2015 Roche Diabetes Care
	1 package insert	results obtained with this test cassette for control solutions can also be traced back to the NIST standard.	Concentrations of galactose in the blood equal to or greater than 20 mg/dL (equal to or greater than		ACCU-CHEK and ACCU-CHEK MOBILE are trademarks of Roche.
	Additional materials required for blood	Detection limit (lowest value displayed): The	1.1 mmol/L) lead to falsely elevated test results.		Roche Diabetes Care GmbH
	glucose testing	detection limit is 10 mg/dL (0.6 mmol/L).	 Do not use when undergoing ceftriaxone treatment. 		Sandhofer Strasse 116
	Accu-Chek Mobile meter with User's Manual Langing daying and langets	Measuring interval: The method is linear within the	Ceftriaxone in the blood may lead to falsely lowered		Germany
	Lancing device and lancets	interval from 10 to 600 mg/dL (0.6–33.3 mmol/L).	test results.		www.accu-chek.com

PDF Approved according to "SOP 04.07.04 - elwo"! The Signature will not be printed! LAY LAN PM LAB TW

8/13/2015 11:56:51 AM

07541538001(01) V1/R1 (Black) - 2015-08 M_Version07_Rev0

pi_07541538001_01_EN(ROW)RU.indd 1

Falzgitter 420x297 / 70x37 Falzschema C TITEL (EN) (RU) **ACCU-CHEK** Mobile 07141254 TITEL (EN) (RU) Мобайл 07141254

(RU) Тест-кассета

Перед тем, как приступить к измерению глюкозы крови с помощью этой тест-кассеты, внимательно ознакомьтесь с информацией, изложенной в этой инструкции-вкладыше и руководстве пользователя глюкометра Accu-Chek Mobile (Акку-Чек Мобайл). В руководстве пользователя содержится вся информация, необходимая для проведения измерения. Если у вас возникли вопросы, обратитесь в Информационный центр.

Данная инструкция-вкладыш содержит 2 следующих символа:

Этот символ указывает на возможность получения травм и причинения вреда вашему здоровью.

этот символ обращает ваше внимание **на важную информацию.**

Назначение

- Тест-кассета предназначена для количественного определения уровня глюкозы свежей капиллярной крови.
- Тест-кассету можно использовать только с глюкометрами Акку-Чек Мобайл и только вне организма пациента
- Система, состоящая из глюкометра и тесткассеты, пригодна только для проведения самоконтроля. Пациенты с диабетом могут контролировать уровень глюкозы крови самостоятельно.
- Систему нельзя использовать для постановки или исключения диагноза диабет.

AKKY-YEK®

- Самоконтроль не заменяет обращения к лечащему врачу. Прежде чем проводить измерение глюкозы крови самостоятельно, обратитесь к лечащему врачу для получения подробных инструкций. Лечащий врач вместе с вами определит индивидуальный, подходящий для вас диапазон уровня глюкозы крови.
 - Храните систему контроля уровня глюкозы крови и все принадлежности вне досягаемости для детей до 3 лет. При проглатывании мелких деталей (например, крышек, насадок или т. п.) существует риск удушья.
- допускается для любых компонентов. находящихся в упаковке. Утилизируйте использованные тест-кассеты согласно требованиям действующего законодательства.
- Получаемые с помощью этой тест-кассеты результаты измерения уровня глюкозы крови соответствуют показателям уровня глюкозы в плазме в соответствии с рекомендациями Международной Федерации по Клинической Химии и Лабораторной Медицине (IFCC) [1]. Таким образом, глюкометр определяет показатели уровня глюкозы крови в плазме, несмотря на то, что на тестовое поле всегда наносится цельная кровь.

Нормальный уровень глюкозы крови не страдающего диабетом взрослого человека не превышает натощак 5,6 ммоль/л (100 мг/дл). Диабет у взрослых диагностируется, если уровень глюкозы крови натощак, подтвержденный двумя измерениями, составляет 7.0 ммоль/л (126 мг/дл) или выше [2, 3, 4]. Состояние взрослых с показателем уровня глюкозы крови натощак от 5,6 до 6,9 ммоль/л (от 100 до 125 мг/дл) определяется как нарушение гликемии натощак (преддиабетическое состояние) [2]. Кроме этих показателей имеются и другие

критерии диагностирования диабета. Обратитесь к лечащему врачу, чтобы определить, страдаете вы лиабетом или нет

Содержимое упаковки

- 1 или 2 тест-кассеты
- 1 инструкция-вкладыш

Тринадлежности, необходимые для еления уровня глюкозы крови

- Глюкометр Акку-Чек Мобайл с руководством пользователя
- Устройство для прокалывания кожи и ланцеты

Для определения уровня глюкозы крови глюкометру требуется 0,3 мкл крови (1 мкл (микролитр) = тысячная миллилитра). Время измерения составляет около 5 секунд (время измерения зависит от концентрации глюкозы крови).

Правильное хранение и применение тест-кассет

- Неправильное хранение или применение тесткассет может привести к получению неверных результатов измерения. Неверные результаты измерения могут привести к неверным рекомендациям по лечению и причинению серьезного вреда здоровью.
- Храните тест-кассеты только в закрытом пластмассовом футляре.
- Храните тест-кассеты при температуре от +2 до +30 °C в сухом месте, защищенном от воздействия прямых солнечных лучей
- При хранении тест-кассеты в холодильнике необходимо подождать, пока закрытый пластмассовый футляр не достигнет температуры окружающей среды. Доставайте тест-кассету из пластмассового футляра, только когда пластмассовый футляр нагреется до температуры окружающей среды. Это позволит предотвратить образование конденсата в тест-кассете.
- Тестовые поля чувствительны к влажности воздуха. При транспортировке тест-кассеты должны находиться в закрытом пластмассовом футляре или в глюкометре. После вскрытия пластмассового футляра тест-кассета должна быть использована в течение 90 дней (срок использования). Использование тест-кассеты по истечении 90 дней может привести к получению неверных результатов измерения уровня глюкозы
- Не используйте тест-кассету, если пластмассовый футляр или защитная пленка повреждены.
- Храните начатую, извлеченную из глюкометра тест-кассету в сухом темном месте.
- Температура, при которой можно проводить измерение, должна составлять от +10 до +40 °C.

Каждое тестовое поле содержит индикаторные реагенты. После нанесения крови на тестовое поле. энзим глюкозодегидрогеназа (Mut. Q-GDH 2 ЕС 1.1.5.2) вступает в реакцию с глюкозой крови. В результате происходящей химической реакции цвет тестового поля меняется. Глюкометр измеряет степень изменения цвета и рассчитывает показатель уровня глюкозы крови.

Контрольное измерение

Чтобы убедиться в корректной работе глюкометра и тест-кассеты, а также в правильном обращении с ними, регулярно проводите контрольные измерения с помощью контрольных растворов Акку-Чек Мобайл. Соблюдайте указания, содержащиеся в инструкции-вкладыше к контрольным растворам, а также в руководстве пользователя глюкометра.

Рабочие характеристики системы Акку-Чек Мобайл

Система Акку-Чек Мобайл отвечает требованиям EN ISO 15197:2013 (в отношении тестовых систем для in vitro диагностики – требования к системам контроля уровня глюкозы крови при сахарном лиабете).

Калибровка и прослеживаемость: В качестве калибровочного средства при калибровке системы (глюкометр и тест-кассета) используется цельная кровь с различными концентрациями глюкозы. Референтные значения определяются гексокиназным методом, который калибрируется методом ID-GCMS. Метод ID-GCMS является методом наивысшего метрологического качества (порядка) и соответствует первичным требованиям стандарта NIST (traceable (прослеживаемость)). Эта цепочка позволяет подтвердить соответствие требованиям стандарта NIST результатов измерения с контрольными растворами, проведенных с помощью тест-кассеты. Нижняя граница показаний (самый низкий отображаемый показатель): Нижняя граница показаний составляет 0,6 ммоль/л (10 мг/дл).

Интервал измерения: Метод является линейным в интервале от 0,6 до 33,3 ммоль/л (от 10 до 600 мг/дл).

Новинка: рабочие характеристики отвечают стандарту EN ISO 15197:2013

Данные о точности системы при концентрации глюкозы ниже 5,55 ммоль/л (ниже 100 мг/дл)

в пределах ± 0,28 ммоль/л	в пределах ± 0,56 ммоль/л	в пределах ± 0,83 ммоль/в		
(в пределах	(в пределах	(в пределах		
± 5 мг/дл)	± 10 мг/дл)	± 15 мг/дл)		
120/156	154/156	155/156		
(76,9 %)	(98,7 %)	(99,4 %)		

Данные о точности системы при концентрации глюкозы равной или выше 5,55 ммоль/л (равной или выше 100 мг/дл)

,			
в пределах ± 5 %	в пределах ± 10 %	в пределах ± 15 %	
216/444 (48,6 %)	410/444 (92,3 %)	442/444 (99,5 %)	

Данные о точности системы при концентрации глюкозы от 1,6 ммоль/л (29,3 мг/дл) до 27,6 ммоль/л (497,0 мг/дл)

в пределах ± 0.83 ммоль/л или в пределах ± 15 % (в пределах \pm 15 мг/дл или в пределах \pm 15 %) 597/600 (99,5 %)

Точность воспроизведения:

	Среднее значение					
	[мг/дл]	43	75	128	214	326
	[ммоль/л]	2,4	4,2	7,1	11,9	18,1
	Стандартное отклонение					
	[мг/дл]	2,0	2,3	3,1	4,3	7,1
	[ммоль/л]	0,1	0,1	0,2	0,2	0,4
	Коэффициент вариации					
	[%]			2,4	2,0	2,2

Промежуточная точность:

	Среднее значение	Среднее значение						
	[мг/дл]	44	126	390				
	[ммоль/л]	2,4	7,0	21,6				
	Стандартное откло	Стандартное отклонение						
	[мг/дл]	2,0	2,8	7,8				
	[ммоль/л]	0,1	0,2	0,4				
	Коэффициент вариации							
	[%]	4,6	2,2	2,0				

Оценка точности пользователем:

Исследование относительно уровня глюкозы в образце, полученном из капиллярной крови, взятой из кончика пальца, проведенного среди 104 обычных людей, показало следующие результаты:

- При концентрации глюкозы ниже 5,55 ммоль/л (ниже 100 мг/дл) 100 % результатов измерений нахолилось в пределах + 0.83 ммоль/л (в пределах ± 15 мг/дл) результатов измерений, полученных лабораторным методом.
- При концентрации глюкозы равной или выше 5,55 ммоль/л (равной или выше 100 мг/дл) 96,8 % результатов измерений находилось в пределах ± 15 % результатов измерений, полученных лабораторным метолом.

Істочники ошибок, которые могут привести к олучению неверных результатов измерения

- Неверные результаты измерения могут привести к неверным рекомендациям по лечению и причинению серьезного вреда злоровью.
- Парентеральное введение галактозы и галактоземия могут привести к получению неверных (завышенных) результатов измерения. Концентрация галактозы в крови равная или выше 1.1 ммоль/л (равная или выше 20 мг/лл) ведет к неверным (завышенным) результатам
- Не применять при лечении цефтриаксоном. Цефтриаксон в крови может привести к неверным (заниженным) результатам измерения.
- При сниженной периферической циркуляции крови капиллярная кровь при определенных обстоятельствах может неверно воспроизводить физиологический показатель уровня глюкозы крови. Это возможно при следующих обстоятельствах: тяжелая дегидратация, обусловленная лиабетическим кетоацилозом гипергликемией/гиперосмолярным состоянием без кетоацидоза, гипотонией, шоком, декомпенсированной серлечной недостаточностью 4-го класса по NYHA или окклюзионным поражением периферических артериальных сосудов
- Можно использовать кровь с показателями гематокрита в диапазоне от 25 до 55 %.

Минимальное содержание на см² на момент изготовления

Варианты мутации хинопротеина глюкозодегидрогеназы (Mut. Q-GDH 2, 6.7 U модифицированный вариант EC 1.1.5.2), acinetobacter spec. 0,40 мкг Пирролохинолин-хинон

Бис-(2-гидроксиэтил)-(4гидроксииминоциклогекса-2,5-8,5 мкг диенилиден) хлорид аммония 2,18-фосфомолибденовая кислота, 88 мкг натриевая соль 0,17 мг Стабилизатор Вещества, не вступающие в реакцию 2,1 мг

Последняя редакция

2015-01

000 «Рош Лиагностика Рус» Бизнес-центр «Неглинная Плаза» **Москва. Россия.** 107031 Информационный центр: 8-800-200-88-99 (звонок бесплатный для всех регионов России)

E-mail: info@accu-chek.ru Адрес в Интернете: www.accu-chek.ru

Список литературы

- [1] D'Orazio et al.: Approved IFCC Recommendation on Reporting Results for Blood Glucose (Abbreviated): Clinical Chemistry 51:9, 1573-1576, 2005
- [2] American Diabetes Association: Standards of Medical Care in Diabetes-2012. Diabetes Care 35 (Suppl. 1), S11-S63, 2012 IDF Clinical Guidelines Task Force, Globa guideline for Type 2 diabetes. Brussels: International Diabetes Federation, 2005
- Definition and diagnosis of diabetes mellitus and intermediate hyperglycemia: report of a WHO/IDF consultation. WHO, Geneva 2006 (ISBN 92 4 159493 4, ISBN 978 92 4 159493 6)



Смотрите инструкцию-вкладыц



указания по технике безопасности, приведенные в инструкции-вкладыше этого Ограничение температуры



(хранить при) Использовать до



Срок хранения тест-кассеты после вскрытия пластмассового футляра 90 дней

Утилизация с бытовыми отходами допускается для любых компонентов

находящихся в упаковке. Утилизируйте использованные тест-кассеты согласно требованиям действующего законодательства.

Производитель

Номер партии

REF

Каталожный номер

LOT IVD

Для in vitro диагностики

GTIN Международный торговый

идентификационный номер

Данный продукт отвечает € 0088 требованиям Европейской Директивы 98/79/ЕС по медицинским устройствам для in vitro диагностики...



ДЛЯ IN VITRO ДИАГНОСТИКИ

© 2015 Roche Diabetes Care

ACCU-CHEK, ACCU-CHEK MOBILE и AKKY-ЧЕК торговые марки фирмы Roche

■ Roche Diabetes Care GmbH

Sandhofer Strasse 116 68305 Mannheim Germany www.accu-chek.com

07541538001(01)

Roche

pi 07541538001 01 EN(ROW)RU.indd 2 8/13/2015 11:56:52 AM